

# Standup support R2D2

The intelligent solution for elderly people and people with various disabilities to use a toilet independently

**Being helped is Good  
but to help yourself is better**

**Aid number for Germany: 33.40.03.1004**

The Standup support R2D2 is an aid for the use of the toilet for everyone who is restricted in his mobility because of his age or physical disability

Thanks to its continuous regulation of height and incline, a safe and joint-protecting standing up and sitting down is guaranteed.

Who doesn't like to maintain his privacy, his independency and his comfort in spite of temporary or permanent physical disabilities?

## **advantages:**

effortless and easy sit down and get up - variable electrical high and inclinable adjustment - easy operation - modular structure - minimum space required, - everywhere and immediately applicable - no assembly expenditure - solid and safety - adjustable on different bowl heights - achievement with and without bidet



## **Vertical Lift**


With folding handles

particularly suitable for wheelchair users and people with muscle weakness in arms.




# Dimensional Drawing


Standup support with fixed arms



Standup support with folding arms



Vertical lift with folding arms



## Standup support with bidet



This combination offers additional a complete personal care

### Technical data

#### Standup support

max.height: 630-860 mm(front-behind)  
 tilting: to 50 degree  
 connection: 230V/50Hz  
 color: white / black  
 drive control: 24VDC  
 load: 150 kg

#### Vertical Lift

max.height: 860 mm  
 connection: 230V/50Hz  
 color: white / black  
 drive control: 24VDC  
 load: 150 kg

### Extras



Control panel



Folding handles



Toilet seat



Dush bidet

### Contact