

VARIANT: 100, 110, 120, 130 & 400

VELA Salsa 100 with a short ergonomic curved seat with rounded front edge for support and stability of the user.



VELA Salsa 100



VELA Salsa 110

VELA Salsa 110 with an ergonomic curved seat with rounded front edge and raised pommel for support and stability of the user. The pommel also prevents the user from sliding down the seat, when tilted.

VELA Salsa 120 with an ergonomic V-shaped front - promoting an active sitting posture and a healthy back.



VELA Salsa 120

*(The footrest is an accessory)*

VELA Salsa is a series of sit-stand chairs with a central brake. The chair promotes the user's independency at different activities, because of a safe, high workplace. VELA Salsa sit-stand chairs promote good ergonomics and are perfect as a second chair throughout a varying work day with sitting and standing tasks at home and at the workplace.

VELA Salsa provides extra safety and stability for the user. A tilted seat supports a healthy posture with an open hip angle and free movement of shoulders and arms e.g. when doing kitchen tasks and working up high.

Stable armrests and a central brake secures safety while getting up and sitting down in the chair. This is an advantage for users with mobility and balance problems. Brake and armrests also provide support when changing the sitting position.

VELA Salsa chairs allows an extra low sitting height, e.g. for the user to help children getting dressed. The chair also comes in a high version, which is optimal when the user is ironing, cooking or for placing dishes in the cupboard.

Adjustment of the seat height, depth and tilt plus backrests and armrests ensures the best conditions for an ergonomic sitting position. This provides support and stability for many activities.

## BENEFITS

- :: VELA Salsa is suitable for sitting as well as standing work, for people with disabilities due to e.g. arthritis, back diagnoses and walking-impairment
- :: Good ergonomics and support ensure an optimum sitting position when working at home and at work
- :: A central brake ensure safety when getting into and out of the chair
- :: It is possible to 'walk' the chair while sitting down

## STANDARD CONFIGURATION

	VELA Salsa 100	VELA Salsa 110	VELA Salsa 120	VELA Salsa 130	VELA Salsa 400
<b>Height adjustment, gas</b>	150 mm (5,9 in), 200 mm (7,9 in) or 260 mm (10,2 in)	150 mm (5,9 in), 200 mm (7,9 in) or 260 mm (10,2 in)	150 mm (5,9 in), 200 mm (7,9 in) or 260 mm (10,2 in)	150 mm (5,9 in), 200 mm (7,9 in) or 260 mm (10,2 in)	150 mm (5,9 in), 200 mm (7,9 in) or 260 mm (10,2 in)
<b>Sitting height</b> The chair can be adjusted in the following intervals: (Tools is necessary)  Note: Measurements in inches are in brackets	<b>min. height at 150 mm:</b> 45-60 cm (17,7-23,6 in) <b>min. height at 200 mm:</b> 50-70 cm (19,7-27,6 in) <b>min. height at 260 mm:</b> 60-86 cm (23,6-33,6 in)	<b>min. height at 150 mm:</b> 45-60 cm (17,7-23,6 in) <b>min. height at 200 mm:</b> 50-70 cm (19,7-27,6 in) <b>min. height at 260 mm:</b> 60-86 cm (23,6-33,6 in)	<b>min. height at 150 mm:</b> 45-60 cm (17,7-23,6 in) <b>min. height at 200 mm:</b> 50-70 cm (19,7-27,6 in) <b>min. height at 260 mm:</b> 60-86 cm (23,6-33,6 in)	<b>min. height at 150 mm:</b> 45-60 cm (17,7-23,6 in) <b>min. height at 200 mm:</b> 50-70 cm (19,7-27,6 in) <b>min. height at 260 mm:</b> 60-86 cm (23,6-33,6 in)	<b>min. height at 150 mm:</b> 45-60 cm (17,7-23,6 in) <b>min. height at 200 mm:</b> 50-70 cm (19,7-27,6 in) <b>min. height at 260 mm:</b> 60-86 cm (23,6-33,6 in)
<b>Seat, width x depth</b>	46 x 35 cm (18,1-13,8 in)	46 x 35 cm (18,1-13,8 in)	40 x 35 cm (15,7-13,8 in)	40 x 35 cm (15,7-13,8 in)	43 x 35 cm (16,9-13,8 in)
<b>Seat tilt, (forwards, backwards)</b>	14° (7° forw., 7° backw.)	14° (7° forw., 7° backw.)	14° (7° forw., 7° backw.)	14° (7° forw., 7° backw.)	14° (2° forw.-16° forw.)
<b>Sitting depth</b>	28-35 cm (11-13,8 in)	28-35 cm (11-13,8 in)	28-35 cm (11-13,8 in)	-	-
<b>Backrest, width x height (Ergo)</b>	38 x 22 cm (15x8,7 in)	38 x 22 cm (15x8,7 in)	38 x 22 cm (15x8,7 in)	-	-
<b>Backrest tilt, (forwards, backwards)</b>	22° (7° forw., 15° backw.)	22° (7° forw., 15° backw.)	22° (7° forw., 15° backw.)	-	-
<b>Height adjustment of backrest</b>	13 cm (5,1 in)	13 cm (5,1 in)	13 cm (5,1 in)	-	-
<b>Armrest padding (PU-armrests)</b>	8 x 26 cm (3,1 x 10,2 in)	8 x 26 cm (3,1 x 10,2 in)	8 x 26 cm (3,1 x 10,2 in)	Optional	-
<b>Armrest height, adjustable</b>	0-24 cm (0-9,4 in)	0-24 cm (0-9,4 in)	0-24 cm (0-9,4 in)	Optional	-
<b>Distance between armrests,</b>	38-48 cm (15-18,9 in)	38-48 cm (15-18,9 in)	38-48 cm (15-18,9 in)	Optional	-
<b>Seat rotation</b>	Optional	Optional	Optional	Optional	No
<b>Frame size</b>	53 x 55 cm (20,1x21,7 in)	53 x 55 cm (20,1x21,7 in)	53 x 55 cm (20,1x21,7 in)	53 x 55 cm (20,1x21,7 in)	53 x 55 cm (20,1x21,7 in)
<b>Chair's weight</b>	17 kg (2st 9,5lb)	17 kg (2st 9,5lb)	17 kg (2st 9,5lb)	14 kg (2st 2,9lb)	14 kg (2st 2,9lb)
<b>User weight</b>	125 kg (19st 9,6lb)	125 kg (19st 9,6lb)	125 kg (19st 9,6lb)	125 kg (19st 9,6lb)	125 kg (19st 9,6lb)

MEASUREMENTS AND WEIGHT

All VELA chairs are CE-marked and comply with applicable standards.

VELA Salsa comes as standard in phoenix fabric, which is durable, breathable and very comfortable. VELA Salsa 110 only comes i Zoom fabric. VELA has its own upholstery workshop that always stocks more than 300 types of foams, textiles and colours. and can fulfill most special requests and customisations.

### ACCESSORIES AND SPECIAL SOLUTIONS

All VELA work chairs can be ordered with a large selection of accessories and special solutions for the chair to be adapted to the needs of individual users. VELA offers individual solutions within padding, support and handles.

A complete list of accessories with item numbers, photos and descriptions is to be found at [www.vela.eu](http://www.vela.eu).

#### VELA SALSA 130 & 400

- :: Stable base with a central brake stabilises the user during activities and when getting in and out of the seat.
- :: VELA Salsa 130 has an ergonomic seat with a V-shaped front - promoting an active sitting posture and a healthy back.
- :: VELA Salsa 400 has a saddle seat which promotes an ergonomically and active sitting posture and ensures a healthy curve in the lower back.



VELA Salsa 130



VELA Salsa 400

Examples of functions for VELA Salsa 130 & 400

VELA Salsa 130 and VELA Salsa 400 does not come with back- and armrest